

Blockthermostat

Quality – Made in Germany



Designed to maintain the temperature of biological materials

Features:

- Temperature adjustable between 27 °C and 45 °C
- Homogeneous temperature distribution
- Heating up time only approx. 40 min.
- Use of different heating blocks

Technical data:

- Dimensions thermostat:
290 mm x 175 mm x 110 mm
(W x D x H)
- Power supply:
100 – 240 V AC, 50/60 Hz, 60 W
- Temperature control between
27 °C and 45 °C
- Accuracy: +/- 0.2 °C
- Temperature distribution: +/- 0.2 °C
- Heating up time from 22 °C to 37 °C:
approx. 40 min.
- Safety class II

REF 10719 (without heating block)

REF 13855 (with heating block: REF 10264)

REF 13896 (with heating block: REF 10277)



Rev. 4_01/2017