

Puncture needles and -sets

Quality – Made in Germany



Labotect ovarian biopsy needles and sets are developed for medical use in the field of assisted reproduction

Product features:

- Needles and sets for follicular puncture
- Facet bevel
- Ergonomic handle
- Grinded or corroded echo marking
- Validated ETO-sterilization
- MEA/LAL tested
- Medical device class IIa

Single lumen puncture sets

with **corroded echo marking for aspiration** only:

REF	Needle length (mm)	Outer diameter of the needle (mm / G)	Plug for tubes with ... ml	Length of tubing from needle to plug (mm)	Length of tubing from plug to aspirator (mm)	Luer connection on plug
13441	300	1.4 / 17	14	1300	1000	yes
13444	300	1.4 / 17	26	1300	1000	no
15138	330	1.0 / 19	14	1300	1000	yes
16783	300	1.4 / 17	14	800	800	yes
16785	330	1.0 / 19	14	800	800	yes

with **corroded echo marking for aspiration & flushing**:

REF	Needle length (mm)	Outer diameter of the needle (mm / G)	Plug for tubes with ... ml	Length of tubing from needle to plug (mm)	Length of tubing from plug to aspirator (mm)	Luer connection on flushing tube
13443	300	1.4 / 17	14	1300	1000	yes

Length of tubing for flushing: 1000

16445	300	1.4 / 17	14	750	30	yes
-------	-----	----------	----	-----	----	-----

Length of tubing for flushing: 750

with **grinded echo marking for aspiration** only:

REF	Needle length (mm)	Outer diameter of the needle (mm / G)	Plug for tubes with ... ml	Length of tubing from needle to plug (mm)	Length of tubing from plug to aspirator (mm)	Luer connection on plug
13438	300	1.4 / 17	14	1300	1000	yes
13440	300	1.4 / 17	26	1300	1000	no
13450	300	1.4 / 17	14	500	1800	no
13451*	300	1.4 / 17	14	1300	1000	no
13452*	300	1.4 / 17	14	500	1000	no
13452L*	300	1.4 / 17	14	1000	1000	no

*with firm tubing

Double lumen puncture sets

with **grinded echo marking** for aspiration & flushing:

REF	Needle length (mm)	Outer diameter of the needle (mm / G)	Plug for tubes with ... ml	Length of tubing from needle to plug (mm)	Length of tubing from plug to aspirator (mm)	Luer connection on flushing tube
13463	330	1.6 / 16	14	390	1000	yes

Length of tubing for flushing: 350 mm


16784	330	1.6 / 16	14	800	800	yes
-------	-----	----------	----	-----	-----	-----

Length of tubing for flushing: 800 mm

Puncture needles

with **grinded echo marking** for aspiration only:

REF	Needle length (mm)	Outer diameter of the needle (mm / G)	Plug for tubes with ... ml	Length of tubing from needle to plug (mm)	Length of tubing from plug to aspirator (mm)	Luer connection
13427	300	1.4 / 17	n.a.	n.a.	n.a.	n.a.



Grinded echo marking



Corroded echo marking